INHALT/INDEX/INDICE

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| | PAGE |
|---|-------------|
| <i>THE TECHNIQUE OF TRANSSPHENOIDAL SURGERY</i> | <i>2</i> |
| <i>REFERENCES.....</i> | <i>19</i> |
| <i>INSTRUMENTS ACCORDING TO PROF. DR. BUCHFELDER.....</i> | <i>20</i> |
| <i>RECOMMENDATION FOR COMPLETE SET</i> | <i>26</i> |

| | PAGINA |
|---|---------------|
| TECNICA DE LA OPERACION TRANSESEFENOIDAL..... | 2 |
| BIBLIOGRAFIA | 19 |
| INSTRUMENTOS SEGUN PROF. DR. BUCHFELDER | 20 |
| RECOMENDACION PARA EL SET COMPLETO | 26 |

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english


TRANSSPHENOIDAL SURGERY

The normal pituitary gland and all tumours that are primarily intrasellar are best approached transspheoidally. Even large intra- and suprasellar pituitary adenomas can be radically resected during such an operation, provided that a wide communication between the intra- and extrasellar tumour portions exists.

*The operation can basically be performed with different types of patient positioning. We currently use the set-up originally described by Cushing where the patient lays in the prone position **(1)** and the surgeon stands behind the head, which is slightly extended **(2)**.*

*Generally, microsurgical and intraoperative imaging techniques are applied. The image intensifier is checked for an orthograde projection **(3)**.*

*The nasal mucosa and the vestibulum oris are infiltrated with a 1% xylocaine solution containing 0,05% adrenaline **(4)**.*

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



english [




We start the operation with a few macroscopical manipulations. Firstly, a small mucosal incision in the vestibulum oris (sublabial) is made to create a mucosal tunnel leading to the sphenoid sinus.

The mucosal incision is carried out about one scalpel width from the plica to the upper lip (5), at a right angle to the raphe (6).

The upper lip is then retracted with a Langenbeck and, by use of the raspatorium, the anterior margin of the cartilagenous nasal septum (7) is exposed for some 10 - 12 mm including the anterior nasal spine (8).

The scalpel and the raspatorium are used to sever the connection of the basal nasal cartilage and mucosa at and, in particular, below the anterior nasal spine (9).

A prominent anterior nasal spine is removed with the diamond drill (10).

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

english [

With the scalpel the cartilage is detached from the anterior nasal spine. The medial nasal mucosa is carefully dissected from the cartilage and the much deeper osseous nasal septum. In this way, a mucosal tunnel is created. It is important to preserve the continuity of the medial nasal mucosa.

The dissection is initially carried out unilaterally and paraseptally, and an attempt is made to completely preserve the cartilagenous nasal septum. Deeper, the osseous septum is encountered **(11)**.

At this stage, the external theatre lights might not provide sufficient illumination of the deep mucosal tunnel. The operative microscope is therefore correctly positioned and used to assist the operation further. Under this microscopic observation, a narrow self-holding dissection speculum with round valves is used to retract the mucosal tunnel **(12,13)**, the cartilage is detached from the osseous nasal septum and portions of this, which would otherwise obstruct a direct approach to the sphenoid sinus, are resected.

At this point, dissection continues via a paraseptal, bilateral route. That is, the medial nasal mucosa is bilaterally exposed, using the keel shaped vomer to provide a safe and reliable midline orientation **(14)**. This technique prevents a deviation to paramedian structures, which would carry the danger of major arterial injury.

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


english [

Just anterior to the vomer, the rounded speculum, that is used for the dissection of the mucosal tunnel, is replaced by a black speculum whose blades are curved outside so that they retract the medial mucosa more laterally (15).

The speculum is not inserted into the sphenoid sinus, but opened anterior to it within the nasal cavity. After opening the sphenoid sinus (16), intrasphenoidal septations and as much as possible of the sphenoid sinus mucosa are resected (17), followed by insertion of a few pads impregnated with hydrogen peroxide to facilitate disinfection and haemostasis.

Then the sellar floor (18) is opened with a microdrill (19) and/or a rongeur (20).

For large tumours which have produced a major expansion of the sellar floor, this is generally very easy. In contrast, small microadenomas often present scenarios where the osseous contour of the sella is hardly enlarged. In this latter situation, a thick sellar floor is likely to be encountered.

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


english [

If the sphenoid sinus is only partially pneumatized, the dorsal and sometimes also the lateral portions of the sella turcica will have to be exposed by use of the diamond drill from the clival region, so that a good overview of the entire intrasellar space is obtained **(21)**.

After incision of the basal dura **(22)**, adenoma tissue can be resected from the intrasellar compartment **(23)**.

The use of the operating microscope allows a reliable differentiation between tumorous tissue, which is softer and whitish in colour, and the normal pituitary gland. In the case of small intrasellar adenomas which measure only a few mm in diameter and which are embedded within the normal pituitary gland, it is necessary to perform multiple sectioning to ensure that the possibly irregularly shaped adenoma is actually completely resected. When applied to microadenomas (tumour diameter < 10 mm), this procedure corrects hormonal oversecretion in some 80 - 90 % of cases.

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

english [

For larger adenomas, neoplastic material in the visible portion of the intrasellar space is initially removed by means of curettes (24) and grasping forceps (25), thereby effecting decompression.

Generally, normal pituitary tissue is displaced into the periphery of the sella, most frequently into the dorsal portion or as a preserved thin yellowish layer just below the sellar diaphragm. Normal pituitary tissue is more yellowish and consistent than tumorous material, and it has a fine vascular structure on its surface. These characteristics help in identification and thus preservation of normal pituitary tissue, and hence of normal postoperative pituitary function. A maximal and radical removal of the sellar floor between both cavernous sinuses facilitates the total resection of large adenomas.

The mostly soft, extrasellar tumour portion spontaneously descends into the space which is created by resection of intrasellar tumour and can thereafter also be removed under direct vision. The descent of the diaphragm into the visible intrasellar space is a good verification of the degree of tumour resection provided there is only one smooth arachnoidal fold (26).

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

english [

However, if the diaphragma sellae descends in multiple folds, there might be residual tumour hidden between these and each folded area needs to be probed separately. If the diaphragma sellae does not spontaneously descend into the intrasellar space after resection of intrasellar tumour, the intracranial pressure should be elevated for a few minutes by either PEEP ventilation or by bilateral compression of the jugular veins.

The lateral portions of the sphenoid sinus and the medial walls of the cavernous sinus can be visualised with a mirror or, if endoscopically assisted microsurgery is carried out, with a 30° or 60° endoscope. This is introduced through the speculum and provides a much better visualisation. In this situation, malleable microinstruments allow well-targeted manipulations even in the lateral aspects of the perisellar space.

After resection of the tumour, the opening in the sellar floor and the diaphragma is covered with an autologous graft, such as derived from the fascia lata of the thigh **(27)**. This is particularly important if there was a CSF leakage intraoperatively or when the arachnoidal layer is very thin and fragile. Human fibrin glue is used as an adhesive **(28)**.

In smaller intrasellar lesions, which are throughout covered by normal pituitary tissue, the application of some kind of collagenous material or fibrin sponge is considered sufficient, provided that there was no intraoperative CSF leakage observed. Resorbable material, such as Gelfoam[®] or Surgicel[®], is inserted into the sphenoid sinus **(29)**.

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

english [

The cartilagenous nasal septum is fixed to the nasal floor by resorbable sutures. Resorbable sutures are also used to close the incision in the vestibulum oris (30).


We recommend insertion for 24 hours of a nasal tamponade (31) extending down to the nasopharynx. This should help to keep the repositioned nasal septum in its position and should prevent epistaxis.

Selective resection of this intra- and suprasellar pituitary adenoma (32) is documented by the delayed postoperative MR investigation. Herein, the complete removal of the tumor and the preservation of the pituitary stalk are visible (33). Since in the early period after surgery, artifacts are anticipated which impede the interpretation of postoperative MR scans, it is recommended to perform postoperative MR studies after 2-3 months following surgery.

A variation, which is particularly suitable for large nostrils and in acromegalic patients, involves a mucosal incision being carried out in the cavum nasi via a perinasal-paraseptal approach. The nostril is initially slightly distended. Then, a mucosal incision, some 3-4 mm below the nasal entrance, is performed throughout the extension of the cartilagenous septum. This mucosal incision should be placed within the nose, as far as possible down from the nostril. Thereafter, by subperichondral dissection in the smooth cleavage plane between the cartilagenous septum and the medial mucosa, a mucosal tunnel is created.

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
english [

As a further variant, there is the option of a direct perinasal approach without any mucosal dissection from the septum at all. For this, the speculum is inserted through one nostril into one nasal cavity and guided under fluoroscopic control towards the sella just anterior to the vomer. The medial mucosa is coagulated and incised in front of the vomer. The sphenoid sinus is opened with the diamond drill and the nasal septum fractured by opening the speculum. The sphenoidotomy is extended with the diamond drill and the rongeur to enable a sufficient overview within the sphenoid sinus. This approach is particularly useful in patients whose noses have been operated on previously, such as following nasal septum corrections by ENT surgeons. It is also very useful when a congenital or acquired nasal septum defect already exists.

During an entirely endoscopic transsphenoidal operation, it is recommended to proceed under vision provided by the endoscope through both nostrils. When we use this procedure, the positioning of the instruments is carefully checked by the image intensifier.

english [

Then, the medial nasal mucosa is coagulated and incised just in front of the vomer and a sphenoidotomy is performed bilaterally. On one side the endoscope, which has been fixed by a reliable device, is introduced. On the other side, microinstruments and suckers can be introduced into the sphenoid sinus. The procedure then follows that described above for the open microsurgical technique.

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english [*Professor Dr. Michael Buchfelder is a senior neurosurgeon at the Department of Neurosurgery, University of Erlangen - Nuremberg in Germany. Following medical studies at the Ludwig-Maximilian-University in Munich he spent his neurosurgical residency in Erlangen with Professor Dr. Rudolf Fahlbusch, with whom he shared a long-term cooperation. The centre in Erlangen is highly modern, well-equipped and has an international reputation, and thus provided an ideal environment for Professor Dr. Buchfelder to pursue his specific interest in technical development of pituitary surgery. In addition, his clinical and scientific interests include the biology, pathophysiology, diagnosis and therapy of pituitary adenomas and other, rare space-occupying lesions in the intra and suprasellar region. In this brochure, the sequence of a typical transspheoidal operation of a pituitary adenoma, as typically performed by the author, is presented together with clear descriptions and illustration of the operative equipment used.*

Professor Dr. Buchfelder enjoys membership of the Deutsche Gesellschaft für Endokrinologie, the Deutsche Gesellschaft für Neurochirurgie, the European Neuroendocrine Association, the Endocrine Society, as well as the American Association and the Congress of Neurological Surgeons.

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Cushing H: **The Weir Mitchell lecture: Surgical experiences with pituitary disorders.** JAMA 63, (1914), pp. 1514 -1525.

Fahlbusch R, Buchfelder M: **The transspheoidal approach to invasive sellar and clival lesions, in: Sekhar LM, Janecka IP (eds) Surgery of cranial base tumors.** Raven Press, New York (1993), pp. 337-349



Griffith HB, Veerapen R: **A direct transnasal approach to the sphenoid sinus.** Technical note. J. Neurosurg 66 (1984), pp. 140-142

Hardy J: **Transspheoidal microsurgery of the normal pathological pituitary.** Clin Neurosurg 16 (1969), pp. 185-217

Jho HD: **Endoscopic pituitary surgery.** Pituitary 2 (1999), pp. 139-154


Laws ER: **Transspheoidal approach to lesions in and around the sella turcica, in: Schmidek HH, Sweet WH (eds) Operative neurosurgical techniques.** Grune & Stratton, New York (1983), pp. 327-341

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





Punch for sellar floor, 45° upwards, 18 cm, with extra thin footplate

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
Image

| No. | Qty. | Art.-Nr. | Description |
|-----|------|-----------|--|
| | 1 | 88.540.15 | Container aluminium 580x280x150 mm, bottom non perforated, lid colour blue |
| | 1 | 88.260.13 | Sterilis.drape green 100x130 cm |
| | 2 | 88.710.26 | Identification label blue, with set description |

Tray 1

| | | | |
|---|---|-----------|---|
| | 1 | 88.590.05 | Wire basket 485x255x50 mm |
| | 1 | 88.260.04 | Sterilisation drape green 40 x 60 cm |
| ① | 1 | 06.107.00 | Scalpel handle no. 7 |
| ② | 1 | 11.308.20 | DEBAKEY atraumatic forceps 2,0/2,5 mm |
| ③ | 1 | 46.154.20 | JANSEN nasal forceps bajonet-shaped, 20cm |
| ④ | 1 | 32.714.16 | WILLIGER periosteal, plain, sharp, 16 cm |
| ⑤ | 1 | 40.442.01 | SPURLING forceps lam. 4x10 mm / 18,0 cm |
| | 1 | 40.422.01 | CUSHING forceps lam. 2x10 mm / 18,0 cm |
| ⑥ | 1 | 41.048.85 | BUCHFELDER dissecting speculum 85 mm, closed valves, for transspenoidal approaches |
| ⑦ | 1 | 41.049.71 | LOYO speculum 70 mm, black |
| | 1 | 41.049.91 | LOYO speculum 90 mm, black |
| | 1 | 41.049.11 | LOYO speculum 110 mm, black |
| ⑧ | 1 | 40.518.02 | Punch for sellar floor, 2 mm, 45° upwards, 18 cm, with extra thin footplate |
| | 1 | 40.518.03 | Punch for sellar floor, 3 mm, 45° upwards, 18 cm, with extra thin footplate |
| ⑨ | 1 | 40.313.23 | FERGUSSON suction tube 2,3 mm/7 Charr. with suction stop |
| | 1 | 40.313.30 | FERGUSSON suction tube 3,0 mm/9 Charr. with suction stop |
| | 1 | 40.313.40 | FERGUSSON suction tube 4,0 mm/12 Charr. with suction stop |
| | 1 | 40.314.23 | FERGUSSON suction tube 2,3 mm/7 Charr. without suction stop |
| | 1 | 40.314.30 | FERGUSSON suction tube 3,0 mm / 9 Charr. without suction stop |
| | 1 | 40.314.40 | FERGUSSON suction tube 4,0 mm/12 Charr. without suction stop |
| ⑩ | 1 | 90.190.30 | Bipolar cord 3 m, flat pin connector for Erbe/ Select/Martin (optionally available for other units) |
| | 1 | 90.603.22 | Bipolar forceps 220 mm, 1,0 mm blunt, bayonet-shaped, angulated up, flat pin connector |
| | 1 | 90.600.22 | Bipolar forceps 220 mm, 0,3 mm pointed bayonet-shaped, flat pin connector |
| | 1 | 90.602.22 | Bipolar forceps 220 mm, 1,0 mm blunt bayonet-shaped, flat pin connector |
| | 1 | 90.602.18 | Bipolar. forceps 180 mm, 1,0 mm blunt bayonet-shaped, flat pin connector |

Image

| No. | Qty. | Art.-Nr. | Description |
|-----|------|-------------|--|
| ⑪ | 2 | 21.029.25 | KOCHER-LANGENBECK retractor 25 x 6 mm |
| ⑫ | 2 | 26.616.12 | HEGENBARTH clip applying forceps, 12,5 cm |
| | 1 | 26.660.16 | MICHEL clips 16 x 3 mm |
| ⑬ | 1 | 08.835.12 | Universal wire scissors TC (to remove clips) |
| ⑭ | 20 | 14.112.09 | TOHOKU towel clamp 10 cm, blunt |
| ⑮ | 1 | 12.221.12 | HALSTED mosquito forceps, curved, 12 cm |
| ⑯ | 2 | 12.310.22 | ROCHESTER PEAN hemostatic forceps, straight, 22 cm (used as a sponge/swab forceps) |
| ⑰ | 2 | 12.311.16 | ROCHESTER PEAN hemostatic forceps, curved, 18 cm (used as a sponge/swab forceps) |
| ⑱ | 2 | 24.364.18 | TOENNIS needleholder TC, 18 cm |
| ⑲ | 1 | 07.197.16TN | LEXER fino scissors blunt/blunt, 16 cm, curved |
| | 1 | 07.283.18TN | METZENBAUM fino scissors, 18 cm, curved |
| ⑳ | 2 | 88.130.06 | Round bowl, stainless steel, Ø 60 mm |

Tray 2

| | | | |
|---|---|-----------|--|
| | 1 | 88.590.03 | Wire basket 485 x 255 x 30 mm |
| | 1 | 88.260.04 | Sterilis.drape green 40 x 60 cm |
| ⑳ | 1 | 41.404.01 | BUCHFELDER diss. rigid, sharp, straight, 26 cm |
| | 1 | 41.404.02 | BUCHFELDER diss. rigid, blunt, straight, 26 cm |
| | 1 | 41.404.05 | BUCHFELDER dura hooklet, bayonet-sh., 26 cm |
| ㉑ | 1 | 41.404.07 | BUCHFELDER pituitary scoop, 26 cm |
| | 1 | 41.404.09 | BUCHFELDER spoon curette sharp, angulated upwards, bayonet-shaped, 26 cm |
| | 1 | 41.404.13 | BUCHFELDER curette blunt 3 mm, downwards, bayonet-shaped, 26 cm |
| | 1 | 41.404.14 | BUCHFELDER curette blunt 5 mm, downwards, bayonet-shaped, 26 cm |
| | 1 | 41.404.18 | BUCHFELDER curette sharp 3 mm, upwards, bayonet-shaped, 26 cm |
| | 1 | 41.404.19 | BUCHFELDER curette sharp 5 mm, upwards, bayonet-shaped, 26 cm |
| | 1 | 41.404.20 | BUCHFELDER curette sharp 8 mm, upwards, bayonet-shaped, 26 cm |
| | 1 | 41.400.52 | LOYO crystal pituitary micromirror, malleable bayonet-shaped, 26 cm |
| ㉒ | 2 | 40.882.16 | NICOLA micro forceps, sharp, 2,5 mm, 16,5 cm |
| ㉓ | 1 | 40.884.16 | NICOLA micro scissors, straight 16,5 cm |
| | 1 | 40.884.15 | NICOLA micro scissors, left 16,5 cm |
| | 1 | 40.884.17 | NICOLA micro scissors, right 16,5 cm |
| ㉔ | 3 | 22.220.15 | Irrigation cannula buttoned, 1,5 x 80 mm LL |
| | 2 | 22.221.12 | Irrigation cannula blunt, 1,2 x 140 mm LL |

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