

6200

Manual Pneumatic Recline with
Power Hydraulic Base Examination Chair

800.735.0357
haag-streit-usa.com



Standard Articulating Headrest

Chair Top and Seat
Crafted From 14-ply Birch

Unique Effortless Recline
Mechanism Supports
Patient Weight

Membrane Switches
with Auto Return

Reclines to the Table Position

All Chair Surfaces Support Patient's Weight

6200

Manual Pneumatic Recline with Power Hydraulic Base Examination Chair

800.735.0357
haag-streit-usa.com

• Ergonomic Features

Reliance chairs are engineered for efficiency and your comfort.

Small footprint allows for increased patient access.

Powered hydraulic base with counterbalanced manual recline.

Dual side controls eliminate the need to reach over a patient to access instrumentation.

Rotation locks allow you to swivel chair 330 degrees from both sides.

Armrests lock in place and release for easy patient access and security.

Headrest adjusts with a single hand and locks into place. Pillow and alternate headrests also available.

Large hinged footrest supports patient's full weight for safe ingress and egress.

Detached footswitch permits hands-free operation of chair base function.

Recline from either side of chair.

• Long-Lasting Membrane Switches

Backlit membrane switches, which last 10-times longer than conventional switches eliminate need for high voltage wiring.

- **Base** raise and lower functions.
- **Auto return switch** brings chair back to lowered position.

• Materials, Finishes and Upholstery

- Long lasting hand-sewn upholstery with high-density foam lasts for years.
- Removable cushions for easy cleaning, or replacement.
- Base available in #20 Pearl.
- 14 standard colors and unlimited custom selections. Visit www.Reliace-Medical.com to review selections.
- Call Reliance at 1-800-735-0357 for a color card.

• The Mark of Quality

The Model 6200 has passed rigorous laboratory testing, which assures quality and safety for you and your patients, achieving CSA, C/US certificate of compliance.

• Tech Specs

Choice of two base heights: High: 21¼" to 33". Low: 19¾" 28¼". Base 22⁵/₁₆ x 24¹/₂. Depth: Upright (without headrest): 40". Table: 60". Seat: 19" (wide) x 18". Armrest span: 19". Shipping weight: 355 lbs. Electrical: 60/50 Hz, 6.3 amperes. Optional electric: 230 volts, 50 Hz, 4.0 amperes achieving CSA, C/US certificate of compliance.

