

BABY SCALE

PRODUCT NUMBER 8320

Medical, health and well-being



- › Foldable to half its length
- › Integrated carrier handle
- › Lightweight and easy to handle
- › Eight-point measuring system guarantees exact weight evaluation

BABY SCALE

PRODUCT NUMBER 8320



Features

The first baby scale with practical collapsible function for parents, midwives and pharmacists can easily reduce its length by half and is impressively light. Furthermore, the integrated carrier handle and successful beauty case design, make it especially manageable and easy to transport. Even with an unsettled baby, the eight-point measuring system guarantees exact weight evaluation, easily legible on a large numeric display. The baby scale contains an environmentally-friendly auto-off function, leading to longer battery life.

Technical Data

- › Electronic weighing system
- › Eight load cells
- › Switchover to lb/oz units
- › Power supply: 2 x 1.5V AA cells
- › Scale switches off automatically for longer battery life

Design

- › Tare function
- › Steady weight display through filtering and hold function
- › Large 20mm LCD digits for easy reading
- › Small space requirement and easy to transport when folded
- › Folded scale held together by magnets
- › Easy and fast opening: Open with carrier handles and remove the base
- › Plastic housing
- › Patent pending

Ordering information

Order Number	Weighing Range	Divisions	Dimensions	Weight
	kg	g	mm	kg
8320.01.001	15	10	375 x 315 x 105 (closed) 740 x 315 x 53 (open)	1,9

Your Soehnle Professional dealer:

www.soehnle-professional.com