

EXAMINATION LIGHT

Halogen Examination
Light for chair models
980, 710, 7000,
6200 and 5200



Rehanee®

Examination Light

Halogen Examination Light for models 980, 710, 7000, 6200, and 5200

Mounts on virtually any Reliance chair for illumination during exams

Application

The Reliance chair mounted examination lamp can be attached to all Reliance Model 980, 710, 7000, 6200, and 5200 chairs.


When specified at the time of chair order, the examination light comes directly wired to Reliance chair models 980, 710, 7000, 6200, and 5200. The light on/off switch is located on the head of the examination light.

Dimensions and Specifications

Gooseneck length: 30". Base housing: 3"x 3" 16"x 2" 16", Power cord length: 96"

Electrical requirements: 115 volts, Approximate shipping weight: 6 pounds

Color: Black



Precision by Tradition