

FX-920

The World's First Full-Power Tilt Exam Chair



800.735.0357
haag-streit-usa.com



Programmable Membrane Switches

Ergonomic Design Saves Space

Power Hydraulic Base Lifts Up To 600 Lbs

Removable Upholstery Aids Cleaning

Adjustable Upholstered Armrests Lock In Place

© 2012 Reliance Medical Products. All Rights Reserved.

FX-920

Patented Full-Power Tilt Exam Chair

800.735.0357
haag-streit-usa.com



● Ergonomic Features

Memory programming cuts wasted movement and improves efficiency.

Dual-side operation controls improve access to patient.

Adjustable headrest and armrests lock with a single hand.

Power tilt eliminates physical recline motion.

Dual-rotation locks allow you to swivel chair 350° without moving.

Large hinged footrest safely supports patient's full weight.

Detached footswitch allows for hands-free patient positioning.

● Power Hydraulic Base and Reclining Top

Finger switches and a corded footswitch raise, lower or tilt chair into position without breaking patient contact, or jolting them. Auto-return button returns patient upright.

● Ultra-Durable Membrane Switches

Backlit membrane switches last 10-times longer than their mechanical counterparts. They're also safer because they run at low voltages. A programmable safety function prevents unwanted operation of the chair.

● Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning or replacement.
- Base available in #20 Pearl.
- Choose from 14 standard vinyl colors and endless custom selections – call Reliance at 1-800-735-0357 for a color card.
- All materials are fire retardant.

● The Mark of Quality

The FX-920 has passed all rigorous Laboratories tests, affirming its safety and durability, achieving CSA, C/US certificate of compliance.

● Tech Specs

Base: 22⁵/₁₆" , width 24¹/₂". Two base heights. High: 24³/₄" lowered and 35³/₈" raised. Low: 22¹/₄" lowered and 30³/₄" raised. Depth upright (excluding headrest): 41"; Reclined: 57". Seat: 20" wide by 15" deep. Distance between armrests: 19 1/2". Shipping weight: 370 lbs. Standard electrical: 110 volts, 60/50 Hz, 6.3 amperes. Optional electrical: 230 volts, 50 Hz, 4.0 amperes.