

1858  
**150**  
2008  
150 Years Haag-Streit

# Goldmann goes digital.

Applanation Tonometer AT 900® D



Goldmann  
Applanation  
Tonometry

**HS** HAAG-STREIT  
INTERNATIONAL

Precision by Tradition

# The digital Goldmann— get your tonometer connected

- True Goldmann applanation tonometer
- Easy to read digital IOP display
- Direct connection to IT system
- Automatic L/R eye recognition
- Proven to have best reproducibility




Goldmann applanation tonometry remains the gold standard for accurate and reproducible IOP measurement.


The digital Goldmann tonometer AT 900 D pairs the positive experience of many thousands of Goldmann tonometer users with the opportunities of state of the art digital technology. Wireless transfer of the IOP readings and left/right eye position allows direct integration of the tonometer into an existing patient information system thus facilitating an ergonomic workflow.

Readings can be used in your existing EMR system or in EyeSuite Tonometry, a Haag-Streit software application included in the optional connection kit.


Goldmann applanation tonometry is simple and reproducible. The Haag-Streit high precision manufacture guarantees continuous reliable performance of the tonometer.



AT 900° D Mod. R



AT 900° D Mod. T



AT 900° D Mod. BQ

Tonometer model			
Slit lamp model	Mod. R	Mod. T	Mod. BQ
BD 900° / BC 900°*	●		
BM 900° / BMV 900°*	●	●	
BP 900°		●	●
BQ 900°		●	●

\* Discontinued slit lamp model

From the original Goldmann manufacturer

## HAAG-STREIT AG

Gartenstadtstrasse 10  
CH-3098 Köniz/Bern, Switzerland  
Tel. 031 978 01 11  
Fax 031 978 02 82  
info@haag-streit.com  
www.haag-streit.com



Precision by Tradition