

5005

Baby height rod

S *Digital and analogue readings*

- › Easy to use, quick and precise
- › With the hold function, the reading remains on the display
- › LCD display, digit height 16 mm
- › Auto-off after 1 minute saves battery power
- › Units in mm and „inch“
- › Incl. adapter for attachment to baby scale 7725 (included in scope of delivery)
- › Battery operation with 2 x 1.5 V AAA batteries (included in scope of delivery)

Order number	5005.10.001
Measuring range	35 – 80.1 cm
Division	1 mm
Dimensions	630 x 290 x 70 mm
Scale weight	60 g

