

980

Premium Full-Power Exam/Procedure Chair

800.735.0357
haag-streit-usa.com



Tufted Removable Upholstery

Armrests Lock in Place

Programmable Safety Function Prevents Unwanted Operation of the Chair

Rounded Base for Easy Access

Footrest Supports Patient's Weight

980

Premium Full-Power Exam/Procedure Chair

800.735.0357
haag-streit-usa.com



● Ergonomic Features

Memory programming cuts wasted movement.

Dual side controls eliminate the need to reach over a patient to access switches.

Rotation locks allow you to swivel chair 330 degrees from both sides.

Armrests lock in place and release for easy patient access and security.

Auto return button returns chair to upright and lowered position.

Large hinged footrest supports patient's full weight for safe ingress and egress.

Detached footswitch allows for hands-free patient positioning.

Headrest adjusts with a single hand and locks into place. Pillow and alternating headrests also available.

● Powered Hydraulic Base and Reclining Top

Keep contact with patient with finger switches and a corded footswitch that raises, lowers or reclines the chair and auto return.

● Long-Lasting Membrane Switches

Backlit membrane switches, which last 10-times longer than conventional switches, contribute to the 980's durability. Eliminate need for high voltage wiring.

- **The programmable safety function** prevents children from operating chair.
- **The memory programmable function** automatically positions chair to preferred exam position.

● Materials, Finishes and Upholstery

- Long lasting hand-sewn upholstery with high-density foam lasts for years.
- Removable cushions for easy cleaning or replacement.
- Selection of 14 standard colors and dozens of custom selections.
- Base available in #20 Pearl.
- Cast iron bases and hinges improve stability and durability.
- Call Reliance at 1-800-735-0357 for a color card.

● The Mark of Quality

The Model 980 has passed rigorous Laboratories testing, which assures quality and safety for you and your patients, achieving CSA, C/US certificate of compliance.

● Tech Specs

Diameter: 24¼", Depth of chair (without headrest) table position: 60", Seat: 18¾" wide, 18¼" deep. Two base heights. HIGH: 23⁵/₈" lowered, 34¼" raised. LOW 21¹/₈" lowered, 29⁵/₈" raised. Distance between armrests: 19" Shipping Weight: 325 lbs. Standard electric: 120 volts, 50/60 Hz, 6.3 amperes total. Optional electric: 230 volts, 50 Hz, 4.0 amperes.