

Original Slit Lamp BM 900®



Original Slit Lamp BM 900®

«When I realised the beauty of what I saw through this slit lamp, I knew that I would become an ophthalmologist.»

This is an authentic statement by Prof. Dr. Georg Eisner, Berne, Switzerland.



The Haag-Streit BM 900® has been the standard in modern Slit Lamp microscopy for more than 40 years – and it still is! Over 60'000 professionals all over the world who own a BM 900® agree, that it was the best investment they ever made for their practice.

Technical Specifications

Irradiance at max. intensity (7.5 V lamp) and max. aperture	0.07 mW/cm ² (300 nm – 400 nm)
	180 mW/cm ² (380 nm – 700 nm)
According to ISO 10939	100 mW/cm ² (700 nm – 1100 nm)

Accessories for the Slit Lamp



Coolpix® Adapter

The Digital-Photo Adapter DA 01 accepts Nikon Coolpix® digital cameras.



Focus attachment

When examining the fundus, the slit may be out of focus and some precise diagnosis becomes difficult. This attachment allows to refocus the slit on the fundus by adjusting the position of the mirror.



Eye-Pieces

A wide range of eye-pieces such as the standard 10x and 16x, wide angle, length and angle measuring eye-pieces are available.

Depth Measuring Attachments

A simple and accurate attachment for measurement of corneal thickness and anterior chamber depth. No I is for measurements up to 1.2 mm. No II for measurements up to 6 mm.



Diffusor

Providing an improved general overview of the eye. Ideal for contact lens fitting because a larger area, beyond 8 mm of diameter, will be illuminated.



Goldmann Applanation Tonometer 900.4.1

Easily mounted on the guide plate for immediate use. When not in use, it is safely stored in the accessory box.



Goldmann Applanation Tonometer 900.4.2.

This instrument remains permanently on the slit lamp and is swung forward in front of the microscope when taking the IOP.

A good slit lamp deserves good contact lenses and HAAG-STREIT has a Diagnostic Contact Lens for every application. Developed in association with Professor Goldmann. Ask your distributor for a complete catalogue.

Here are just three of the HS Contact Lens range:



Retina 145



Goldmann 3-mirror 903



Fundus 901

Motorised table, manual table or diagnostic unit – HAAG-STREIT quality solutions for everyone.



HSM 801

HSM 901

SWING Diagnostic Unit