

# 3000

Manual Pneumatic Exam and Treatment Chair



800.735.0357  
haag-streit-usa.com



# 3000

## Pump Hydraulic Lift and Manual Recline Chair

800.735.0357  
haag-streit-usa.com

### ● Ergonomic Features

**Adjustable padded headrest** locks in place with a single hand.

**Foot bar and armrests** safely support patient's full weight.

**Armrests** rotate to allow patients in or out from both sides of the chair.

**Rotation lock** allows chair to swivel 360° for improved patient access.

### ● Pump Hydraulic Base for Effortless Raise

Hydraulic mechanism allows you to easily raise the patient by pumping the 3000's foot pedal. Automatically lowers when pedal is moved into the down position.

### ● Manual Reclining Top

Chair back reclines up to 90° with a single hand. Convenient release mechanism on both sides of the chair.

### ● Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning or replacement.
- Interchangeable seat and back cushions.
- Choose from 14 standard vinyl colors, or endless custom selections – call Reliance at 1-800-735-0357 for a color card.

### ● Tech Specs

Top and base shipped together for easy installation. Base height: 28" raised, 19" lowered.

Base dimensions: 23" diameter. Overall width: 25". Depth upright: 32½". Depth reclined: 42".

Seat dimensions: 19" by 19". Armrest distance: 20". Shipping weight: 220 lbs.

