

522

Stylish Exam Chair with Small Footprint



800.735.0357
haag-streit-usa.com



522

Stylish Exam Chair Conveniently Slides into Any Lane

800.735.0357
haag-streit-usa.com



Overview

With a contemporary design, the 522 is a perfect fit for a new practice, satellite office, or just about anywhere.

Ergonomic Features

Large padded headrest keeps patients comfortable during exams.

Dual-rotation locks allow you to rotate chair 360° from a single position.

Upholstered armrests rotate up to allow patients to move in and out from either chair side.

Large hinged footrest safely supports patient's full weight, and folds away to save space.

Stationary back allows the 522 chair to fit into any size exam room.

Power Hydraulic Base Boosts Efficiency

The 522 chair's corded footswitch allows you to raise and lower the patient while maintaining contact.

Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Choose from 14 standard vinyl colors and 235 custom selections. Call Reliance at 1-800-735-0357 for a color card.
- Removable cushions for easy cleaning, or replacement.
- Base available in #20 Pearl.
- All materials are fire retardant.

CSA Mark of Quality

The 522 has passed all rigorous CSA tests, affirming its safety and durability. Made in America.

Tech Specs

Low base height: 19½" lowered, 28½" raised. Base: 22⅞" by 24¼". Depth upright: 43½". Seat: 21½" wide by 15" deep. Armrest distance: 21⅝". Shipping weight: 310 lbs. Standard electrical: 120 volts, 60 Hz, 3.0 amperes. Optional: 230 volts, 50 Hz, 3.0 amperes.