

710



Minor procedures  
chair/table with  
**pillow headrest**

**Reliance®**

# 710

Features foldaway headrest and reclines 180° for minor procedures.



## Surgical Headrest

A non-articulating pillow headrest is standard. An optional special articulating pillow headrest features a sliding, removable magnetic pad and a fold-away design for otology. Available on Model 710 only.

## Powered Base And Articulating Top

Use finger switches and corded foot switch to raise, lower or recline chair up to 180° without breaking contact with patient. Eliminates wasted movements. An auto return button automatically returns the chair to the low, upright position.

## Memory Positioning/Safety Switch

A memory switch allows you to instantly pre-select the chair's most used position. The patented programmable safety switch disables side and foot switches which prevents children from operating the chair.

## Membrane Switches

Membrane switches enhance the Model 710's sleek design and contribute to its durability and safety. They last 10-times longer than mechanical switches, eliminate the need for high voltage, and allow for limited programming of chair functions. Back lighting eases work in a dimmed room. Mechanical switches available as option.

## Dual Rotation Locks

Easily accessible levers located on both sides of seat enable practitioner to swivel chair 330° without changing position.

## Large Hinged Footstep

The 8" x 13" footstep supports patient's legs, and provides safe ingress and egress. Tilts out of way. Upholstered pad supports patient's feet for standard procedures.

## Armrests

Armrests lock in place. Release mechanism allows convenient ingress and egress from either side of chair. Padded underneath for comfort.

## Upholstery

Match office or treatment room décor with our selection of 14 easy-to-maintain vinyl colors. Other colors can be supplied at additional cost. Available with either smooth or the designer-style cushions shown on cover. Base is available in two standard enamel colors: #17 Light Beige and #15 Silver Smoke. All materials are fire retardant. To request a color selection card, contact your dealer or call Reliance at 1-800-735-0357.

## The UL Mark Of Quality

The Model 710 chair has passed rigorous testing by Underwriters Laboratories, which assures its quality and safety for you and your patients.

## Dimensions and Specifications

Base dimensions: 22 1/8" x 24" diameter. Overall width: 24 1/4". Overall depth in chair position (excluding headrest): 40". Table position: 64". Seat: 18 1/2" (wide) x 18" (deep). Distance between armrests: 19". Choice of two base heights: HIGH: 22 1/2" lowered, 34" raised; LOW: 20" lowered, 29" raised. Shipping weight: 290 lbs. Electrical: 115 volts, 60 Hz, 8.0 amperes.



Precision by Tradition