

7720



Ophthalmic  
**instrument stand**  
with textured finish

**Rehance**<sup>®</sup>

# 7720

Traditional design. Excellent value.

## The Perfect Match for all Reliance Chairs

The Model 7720 IC (with console) and the 7720 NC (without console) are designed for maximum access to patients and compatibility with all Reliance examination chairs. They include a tray area, a heavy duty scuff pad on the front of the base casting, and chair controls.

## Textured Finish

Tough, textured epoxy finishes are offered in two standard colors. To request a color selection card, contact your dealer or call Reliance at 1-800-735-0357.

## Counterbalanced Slit Lamp Arm

A locking mechanism under the slit lamp arm lets you adjust its height with a single hand. It has a 12" range of travel, and it can be counterbalanced for instruments weighing between 12 and 60 lbs. A switch to raise and lower the chair is available on the slit lamp arm.

## Upper (Refractor) Arm

The refractor arm has a 10" range of vertical travel and accommodates instruments weighing up to 20 lbs. It can be positioned on the support column at the desired height. A third instrument arm with a 20 to 50 lb. capacity is offered as an option. Its vertical travel is 10".

## Instrument Wells (IC only)

The Model 7720-IC is equipped with three instrument wells which hold rechargeable Welch Allyn instruments. A hanger can be provided for an indirect ophthalmoscope.

## Membrane Switches

Backlit membrane switches on the instrument console control the chart projector and the raise, lower and auto down functions of a Reliance examination chair. Membrane switches last 10-times longer than mechanical switches, and increase safety by eliminating high voltages. They also control the overhead lamp, and the binding posts typically used for the indirect ophthalmoscope (IC only).

## Halogen Overhead Lamp

A standard overhead lamp arm on the Model 7720 is equipped with a long-lasting halogen bulb. It provides a bright, direct light source for minor procedures. Low voltage is safer, and the arm adjusts without the need for tools.

## Dimensions and Specifications

Height to top of console: 36". Height to top of pole: 71". Base: 21" x 21". Shipping weight: 292 lbs. Electrical: (Standard; UL Listed) 115 volts, 60Hz, 3 amperes; (Optional) 230 Volts, 50 Hz, 1.5 amperes.

## Options

**Wheel Chair Access:** An optional extended refractor arm and a modified lower slit lamp arm make it possible to deliver instruments to patients in any model of wheelchair. When ordered with this option, the cost may be applied toward the ADA tax credit.

**Eclipse™ Room Light Control:** An infrared system controls all room lighting through instrument moves or console switches. Available as an option on all Reliance low voltage stands. Call authorized Reliance dealer for more information.

**The Model 7725:** For practitioners who only require a slit lamp arm on an instrument stand. This same sturdy construction and features as the 7720 stand without a refractor arm, lamp and pole. Available with (IC) or without (NC) an instrument console.



Model 7720 IC  
With Instrument  
Console



Model 7720 NC  
Without Instrument  
Console

**HS** RELIANCE  
INTERNATIONAL

Precision by Tradition