

# Easy use, optimal imaging.

Imaging Module IM 900



Digital Slit  
Lamp Imaging

 **HAAG-STREIT  
INTERNATIONAL**  
Precision by Tradition


# The intelligent choice for optimal imaging.

- Easy intuitive imaging
- Still imaging and high resolution video
- Freeze Technology – high speed capture and data transfer
- Easy Touch – precise image control
- DFC – depth of field control
- History Trigger – never miss the perfect image
- Compact

Imaging Module IM 900 the fully integrated compact imaging solution for the BQ 900 slit lamp.


Combination of the IM 900™ with EyeCap™ or EyeSuite™ Imaging completes the perfect imaging system.

## FREEZE TECHNOLOGY



**Freeze technology:** Captures at the precise moment you press the trigger.

## HISTORY TRIGGER



**Rewind time and freeze the moment:** The history trigger not only captures one image but records the last few seconds in real time and lets you return to select the optimum image.

## EASY TOUCH



**Control the image:** Easy touch more than just a trigger, allowing simple, ergonomic management of the camera settings.

## DEPTH OF FIELD CONTROL (DFC)



**True aperture:** A unique feature of the Haag-Streit imaging systems, facilitating complete control of the depth of field.



From the leader in slit lamp imaging

**HAAG-STREIT AG**  
Gartenstadtstrasse 10  
CH-3098 Köniz  
Tel. +41 31 978 01 11  
Fax +41 31 978 02 82  
info@haag-streit.ch  
www.haag-streit.com

**HS HAAG-STREIT  
INTERNATIONAL**  
Precision by Tradition