



ClearChart[®]2

Digital Acuity System

Smart, easy acuity.



Reichert's **ClearChart®2** Digital Acuity System builds on the immensely successful ClearChart system. This is a comprehensive solution for your entire acuity chart needs. The simple to use, all-in-one design is sleek and attractive eliminating the need for a PC or additional components.

Featuring the all-new *19 inch*, high-resolution LCD flat panel display, and weighing only 11.25 pounds, the ClearChart 2 is easily mounted to any wall. The custom designed LINUX system is free of fans and hard disk drives, allowing for silent, dependable operation in the exam room. ClearChart 2's convenient, infrared remote control enables effortless access of all tests from virtually anywhere in the room.



Snellen Chart



Contrast Sensitivity



Pediatric Fixation



Patient Education

ClearChart 2 features:

ClearChart 2 incorporates all the features and functionality you need for effective and efficient visual acuity testing.

- Large range of optotypes and tests
- Random sequencing
- Contrast acuity and sine-wave grating contrast sensitivity testing
- Cartoon loop with sound for pediatric fixation
- Adjustable Red/Green for compatibility with any refraction system
- Configure for direct or mirrored viewing.
- Complies with ETDRS protocol
- Patient education slides



ClearChart 2 enhances the refraction experience for both the doctor and the patient. The bright, high-contrast display enables procedures to be performed with the lights on or off. The easy to use remote allows the examiner to switch exam screens, exam modes, and access slides to help educate the patient, by merely pressing a button.

Reichert quality:

The Reichert and ClearChart names ensure the highest quality and reliable instrument that thousands have come to depend on for their refraction needs.

Catalog Number

13760 ClearChart 2 Digital Acuity System with Wall Mount

Specifications

Display	19 inches (48.26 cm) LCD
Max Illumination	220 cd/m ²
Photopic Illumination	85 cd/m ²
Mesopic Illumination	3 cd/m ²
Height	16.91 inches (42.95 cm)
Width	19.88 inches (50.49 cm)
Depth	2.52 inches (6.4 cm)
Weight (unpacked)	11.25 lbs. (5.1 Kg)
Power Requirements	100-240 volts AC 50/60 Hz 93-109 VA
Fuses	T3.15AL250
Refraction Range	6 ft (1.83 m) - 31 ft (9.5 m)

Proudly Made in the USA.

Visit Reichert Technologies at www.reichert.com or contact your Authorized Reichert Technologies Distributor.

Important Notice

For your own protection, instruments should be purchased from and serviced by Reichert Technologies or an Authorized Reichert Technologies Distributor only.

CC2-1010 ©2010 Reichert, Inc. Specifications subject to change without notice. Reichert and ClearChart are registered trademarks of Reichert, Inc. Reichert Technologies, and "Advancing Eye Care. Preserving Sight" are trademarks of Reichert, Inc.

Reichert
TECHNOLOGIES
Advancing Eye Care. Preserving Sight.™

Reichert, Inc.
3362 Walden Avenue, Depew, NY 14043
Toll Free: 888-849-8955 • Tel: (716) 686-4500
Fax: (716) 686-4545 • www.reichert.com

ISO 9001 Certified
Reichert products are designed and manufactured under quality processes meeting ISO 9001 requirements.