

# 7000

## Dual-Purpose Exam Chair and Minor Procedure Table

800.735.0357  
haag-streit-usa.com



# 7000

## Dual-Purpose Exam Chair and Minor Procedure Table

800.735.0357  
haag-streit-usa.com



### ● Ergonomic Features

**Dual-side operation** allows you to control chair from both sides while maintaining patient contact.

**Adjustable concave headrest** provides superior head and neck support, and locks with a single hand.

**Upholstered armrests** lock in place and release on both chair sides for security and comfort.

**Large hinged footrest** supports patient's full weight and tilts out of way to save space.

**Corded footswitch** provides hands-free operation of chair base and top. Dual rotation locks allow you to swivel chair 330°.

### ● Perfect for Minor Procedures

The 7000 chair reclines nearly flat for minor procedures. Non-articulating pillow headrest and articulating surgical headrest are optional.

### ● Power Hydraulic Square Base

Square pedestal base provides crucial patient stability and superior accesability. Finger switches and a corded footswitch raise, lower or recline chair into any position without breaking patient contact. High-and low base models available.

### ● Ultra-Durable Membrane Switches

Programmable backlit membrane switches last 10-times longer than their mechanical counterparts. They're also safer because they run at low voltages.

### ● Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning or replacement.
- Designer-style cushions standard.
- Base available in #20 Pearl.
- Choose from 14 standard vinyl colors and endless custom selections – call Reliance at 1-800-735-0357 for a color card.

### ● CSA Mark of Quality

The 7000 has passed all rigorous testing by CSA certification to assure its safety and durability, achieving CSA, C/US certificate of compliance.

### ● Tech Specs

Base: 22<sup>5</sup>/<sub>16</sub>" x 24<sup>1</sup>/<sub>2</sub>". Overall width: 24<sup>1</sup>/<sub>2</sub>". Overall depth (excluding headrest): 36". Table position: 60". Seat: 19" wide by 18<sup>1</sup>/<sub>2</sub>" deep. Armrest distance: 19". Choice of two base heights. High: 22<sup>3</sup>/<sub>4</sub>" lowered, 33<sup>3</sup>/<sub>8</sub>" raised. Low: 20<sup>1</sup>/<sub>4</sub>" lowered, 28<sup>3</sup>/<sub>4</sub>" raised. Shipping weight: 325 lbs. Standard electrical: 115 volts, 60/50 Hz, 6.3 amperes total. Optional: 230 volts, 50 Hz, 4.0 amperes.